

General Information

Test-Box contains one pair each of **OHROPAX® Classic**, **OHROPAX® Soft** and **OHROPAX® Silicon PINK** earplugs.

OHROPAX® Classic

Earplugs against harmful noise for single use only, to be moulded before use (EN 352-2:2002, EN 2016 / 425).

Notified body for conformity assessment and production control: IFA Institute for Occupational Safety and Health of the German Social Accident Insurance (0121), Alte Heerstraße 111, 53757 St. Augustin, Germany.

Ingredients: Vaseline, paraffin wax, cotton wool

Colour: Blush-pink

Size: Not applicable

Characteristics: Mouldable, skin-friendly wax material. Adapts to any ear without pressure.

Instructions: ON-Ear plug (to be placed in front of the opening of the ear canal)



Before use, completely remove cotton wool.
Use entire earplug – do not split!



Soften the plug by vigorous kneading, form it into a ball (no elongated shape!)



Seal the opening of the ear canal with gentle pressure. **Do not insert into the ear canal!**


Removal: Use your fingers to remove the plug carefully in one piece. Do not use any auxiliary means! The plug can be loosened by slightly pushing the back of the ear.

Sound Attenuation (EN ISO 4869-2)

Frequency (Hz)	63	125	250	500	1000	2000	4000	8000
Mean value (dB)	21,4	20,8	18,8	19,3	25,7	35,0	43,2	38,8
Standard deviation (dB)	6,3	6,2	7,3	6,8	5,1	4,6	7,6	11,6
Assumed protection value APV (dB)	15,1	14,6	11,5	12,5	20,6	30,4	35,6	27,2

SNR* = 22 dB H* = 27 dB M* = 18 dB L* = 14 dB

* Value



Further Information

- Insert the earplugs as described and check that they fit correctly.
- Ensure hygienic use. Your hands, ears and earplugs should be clean and dry!
- The product cannot be cleaned.
- Replace soiled or damaged earplugs or such plugs that have become less effective.
- Wear earplugs at all times in noisy environments.
- Sudden or fast removal of the earplug out of the ear may damage the eardrum.
- This product may be adversely affected by certain chemical substances. Further information can be sought from the manufacturer.
- Warning: Failure to comply with these instructions may significantly impair the protective effect of the product.
- Keep away from small children (danger of swallowing).
- Can be used by children aged six and older; use only under adult supervision.
- Storage: In the box, dry and below 30°C.

Application Areas

For private use: OHROPAX® Classic is ideal for sleep and relaxation due to its high wearing comfort, but also for concentrated work, DIY and loud music. Also keeps moisture and wind away from your ear.

At work:

- OHROPAX® Classic is particularly suitable in the so-called non-traditional noise areas such as open-plan offices and computer centres (Noise level <75 dB (A)).
- OHROPAX® Classic is immediately ready for use. There is no need to wait for expansion times as with foam plugs.
- No sweating and no pressure on the auricle as with capsule hearing protectors.
- No wearing problems for spectacle wearers.
- Can be combined with safety helmets, goggles and capsule hearing protectors.
- Speech and warning signals are still perceived. Sound is attenuated, but not completely excluded.
- OHROPAX® Classic is not suitable for warm, dirty and dust-laden workplaces due to its material.

Suitability

Before using the product in noisy environments, first check whether the plug is suitable for this load. Take into account the special notes, sound attenuation values and please see the [information on sound attenuation values and protective effect](#).

Tolerability

OHROPAX® Classic is considered to be extremely well tolerated. Components with known allergenic potential are not included. Sometimes you have to get used to hearing protection. We recommend increasing the wearing time gradually. In rare cases, you can initially hear your own pulse. This effect usually disappears after a short period of habituation.

General Information

Reusable silicone earplugs against harmful noise, to be moulded before use (EN 352-2:2002, EN 2016/425).

Notified body for conformity assessment and production control: IFA Institute for Occupational Safety and Health of the German Social Accident Insurance (0121), Alte Heerstraße 111, 53757 St. Augustin, Germany.

Material: Silicone
Colour: Pink
Size: Not applicable

Characteristics: Mouldable, pure silicon. The skin-friendly earplugs adapt to any ear without pressure.

Instructions: ON-Ear plug (to be placed in front of the opening of the ear canal)



Soften the plug by vigorous kneading, form it into a ball (No elongated shape!).

Use entire earplug – do not split!



Place the ball in front of the opening of the ear canal and seal it with gentle pressure.

Do not insert into the ear canal.




Removal: Use your fingers to remove the plug in one piece. Don not use any auxiliary means. The plug can be loosened by slightly pushing the back of the ear. If you have long fingernails, be careful not to damage the earplug.

Sound Attenuation (EN ISO 4869-2)

Frequency (Hz)	63	125	250	500	1000	2000	4000	8000
Mean value (dB)	23,5	21,4	19,1	20,3	23,6	30,7	36,4	35,0
Standard deviation (dB)	4,7	3,5	2,8	4,3	4,5	4,7	3,1	6,2
Assumed protection value APV (dB)	18,8	17,9	16,3	16,0	19,1	26,0	33,3	28,8

SNR* = 23 dB H* = 26 dB M* = 19 dB L* = 18 dB

* Value



Further Information

- Insert the earplugs as described and check that they fit correctly.
- Ensure hygienic use. Your hands, ears and earplugs should be clean and dry!
- The earplugs should be regularly inspected for usability as well as for the necessity of cleaning.

- Replace soiled or damaged earplugs or such plugs that have become less effective.
- Cleaning: Light impurities can be wiped off using a cloth and warm water. Wait until completely air-dried.
- Wear earplugs at all times in noisy environments.
- Sudden or fast removal of the earplug out of the ear canal may damage the eardrum.
- **For water protection:** Do not use for diving at depths greater than 1m. Stronger movements when swimming or jumping into the water can cause the ear protection to loosen and get lost. If the ear is inflamed, the eardrum is injured or if tympanic tubes are used, consult your doctor first.
- This product may be adversely affected by certain chemical substances (e.g. hand cream). Further information can be sought from the manufacturer.
- Make sure the plug is large enough not to accidentally enter the ear canal. If removal with the fingers is no longer possible, please consult a doctor (preferably an ENT specialist).
- Warning: Failure to comply with these instructions may significantly impair the protective effect of the product!
- Keep away from small children (danger of swallowing).
- Can be used by children aged six and older; use only under adult supervision.
- Storage: In the box, dry and below 30°C.

Application Areas

For private use: OHROPAX® Silicon PINK is ideal for sleep and relaxation due to its high wearing comfort, but also for concentrated work, DIY and against loud music. It also protects the ear against the entry of water when bathing.

At work:

- OHROPAX® Silicon PINK is suitable for many workplaces, especially in the so-called non-traditional noise areas such as open-plan offices and computer centres (Noise level <75 dB (A)).
- OHROPAX® Silicon PINK is immediately ready for use. There is no need to wait for expansion times as with foam plugs.
- Speech and warning signals are still perceived. Sound is attenuated, but not completely excluded.
- OHROPAX® Silicon PINK is not suitable for dirty and dust-laden workplaces due to its material.

Suitability

Before using the product in noisy environments, first check whether the plug is suitable for this load. Take into account the special notes, the sound attenuation values and please also see the [information on sound attenuation values and protective effect](#).

Tolerability

OHROPAX® Silicon PINK is considered to be extremely well tolerated. Components with known allergenic potential are not included. Sometimes you have to get used to hearing protection. We recommend increasing the wearing time gradually. In rare cases, you can initially hear your own pulse. This effect usually disappears after a short period of habituation.

General Information

Reusable foam earplugs against harmful noise (EN 352-2:2002, EN 2016/425).

Notified body for conformity assessment and production control: IFA Institute for Occupational Safety and Health of the German Social Accident Insurance (0121), Alte Heerstraße 111, 53757 St. Augustin, Germany.

Material: PU foam

Colour: Blush-pink

Size: Ø 6-12 mm

Characteristics: The soft, anatomically shaped earplugs gently expand in the ear canal and offer reliable protection against harmful and disruptive noise. They provide a good fit due to their open-pored, slightly rough surface.

Instructions: IN-Ear plug (to be inserted into the ear canal)



Roll the earplug into a tight cylinder.



Reach over your head with one hand and pull the ear slightly upward with one hand, then insert the rolled earplug deep into the ear canal with the other.



Hold the earplug in place until it has expanded and fits tightly.

Removal: Remove the plug carefully with your fingers.

Sound Attenuation (EN ISO 4869-2)

Frequency (Hz)	63	125	250	500	1000	2000	4000	8000
Mean value (dB)	30,6	32,9	32,5	34,8	34,6	32,5	39,3	41,2
Standard deviation (dB)	6,7	4,9	4,6	6,1	6,5	3,8	3,2	3,7
Assumed protection value APV (dB)	23,9	28,0	27,9	28,7	28,1	28,7	36,1	37,5


SNR* = 32 dB

H* = 31 dB

M* = 29 dB

L* = 28 dB

* Value



Further Information

- Insert the earplugs as described and check that they fit correctly.
- Ensure hygienic use. Your hands, ears and earplugs should be clean and dry!
- The earplugs should be regularly inspected for usability as well as for the necessity of cleaning.
- Replace soiled or damaged earplugs or such plugs that have become less effective.
- For cleaning (up to two times) wipe the surface of the plug (without pressure) with a mild soap solution. Then rinse and wait until completely air-dried. No wringing!
- Wear earplugs at all times in noisy environments.
- Sudden or fast removal of the earplug out of the ear canal may damage the eardrum.
- This product may be adversely affected by certain chemical substances. Further information can be sought from the manufacturer.
- Warning: Failure to comply with these instructions may significantly impair the protective effect of the product!
- Keep away from small children (danger of swallowing).
- Can be used by children aged six and older; use only under adult supervision.
- Storage: In the box, dry and protected from heat.

Application Areas

For private use: If you are exposed to noise in your free time or just want to keep your peace, OHROPAX® Soft are the ideal earplugs.

At work:

- OHROPAX® Soft is suitable for many workplaces, especially in the so-called non-traditional noise areas such as open-plan offices and computer centres (Noise level <75 dB (A)).
- No sweating and no pressure on the auricle as with capsule hearing protectors.
- No wearing problems for spectacle wearers.
- Can be combined with safety helmets, goggles and capsule hearing protectors.
- Speech and warning signals are still perceived. Sound is attenuated, but not completely excluded.

Suitability

Before using the product in noisy environments, first check whether the plug is suitable for this load. Take into account the special notes, the sound attenuation values and please also see the [information on sound attenuation values and protective effect](#).

Tolerability

OHROPAX® Soft is considered to be extremely well tolerated. Components with known allergenic potential are not included. Sometimes you have to get used to hearing protection. We recommend increasing the wearing time gradually. In rare cases, you can initially hear your own pulse. This effect usually disappears after a short period of habituation.